

SAFETY RECALL

Bed Handles, Inc.

**Bedside Assistant[®] Travel Handles[™] and
Adjustable Bedside Assistant[®]
Adult Bed Handles**

When attached to an adult bed, without the use of safety retention straps, the handle can shift out of place, creating a dangerous gap between the handle and the side of the mattress posing a risk of entrapment, strangulation and death.



**Bedside Assistant
Model# BA10W**

**Travel Handle
Model# BA11W**

**Adjustable Bedside Assistant
Model# AJ1**

Sold nationwide through 2007

Immediately stop using the bed handles and contact the firm to receive free safety retention straps to secure the bed handle to the bed frame, along with installation instructions and warning label.

1-800-725-6903

8:30 a.m. to 4:30 p.m. CT Monday – Friday or online

www.BedHandles.com